

## NOBEL TORQUE KIT


Designation: Dispensary

Sterilization Cycle: 4/40 @ 134°C

Wrap Method: Own Container / Steri - Peel

Quantity: 20

1. Nobel Torque Wrench Prosthetic 29165
2. Nobel Torque Wrench Adaptor 29167
3. Nobel Multi-Unit Screw Driver 21mm 29158
4. Nobel Locator Latch Wrench Driver 23mm 08914
5. Nobel Screwdriver Machine UniGrip 30mm 29153
6. Nobel Screwdriver Machine UniGrip 25mm 29152
7. Nobel Screwdriver Machine UniGrip 20mm 29151



### Cleaning:

Washer/Disinfector (instruments cycle)